



# Anti-GGNBP2 (aa 223-320) polyclonal antibody (DPABH-18990)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                                                                                          |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antigen Description</b> | May be involved in spermatogenesis.                                                                                                                                                                      |
| <b>Immunogen</b>           | antigen sequence: VLRAYNILIG ELDCSKEKGY CAALY EGLRC C PHERHIHVC CETDFIAHLL GRAEPEFAGG RRERHAKTID IAQEEVLTCL GIHL YERLHR IWQKLRAE, corresponding to internal sequence amino acids 223-320 of Human GGNBP2 |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                      |
| <b>Source/Host</b>         | Rabbit                                                                                                                                                                                                   |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                                                    |
| <b>Purification</b>        | Immunogen affinity purified                                                                                                                                                                              |
| <b>Conjugate</b>           | Unconjugated                                                                                                                                                                                             |
| <b>Applications</b>        | IHC-P                                                                                                                                                                                                    |
| <b>Format</b>              | Liquid                                                                                                                                                                                                   |
| <b>Size</b>                | 100 µl                                                                                                                                                                                                   |
| <b>Buffer</b>              | pH: 7.20; Constituents: 59% PBS, 40% Glycerol                                                                                                                                                            |
| <b>Preservative</b>        | 0.02% Sodium Azide                                                                                                                                                                                       |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.                                                                                                           |

## GENE INFORMATION

|                        |                                                                                                                                                                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>       | <a href="#">GGNBP2 gametogenetin binding protein 3 [ Homo sapiens ]</a>                                                                                                                                                                                 |
| <b>Official Symbol</b> | GGNBP2                                                                                                                                                                                                                                                  |
| <b>Synonyms</b>        | GGNBP2; gametogenetin binding protein 2; DIF3; LZK1; DIF-3; LCRG1; ZFP403; ZNF403; gametogenetin-binding protein 2; zinc finger protein 403; C3HC4-type zinc finger protein; laryngeal carcinoma related gene 1; laryngeal carcinoma-related protein 1; |
| <b>Entrez Gene ID</b>  | <a href="#">79893</a>                                                                                                                                                                                                                                   |
| <b>Protein Refseq</b>  | <a href="#">NP_079111.1</a>                                                                                                                                                                                                                             |
| <b>UniProt ID</b>      | <a href="#">Q9H3C7</a>                                                                                                                                                                                                                                  |