



Anti-NYNRIN (aa 1476-1604) polyclonal antibody (DPABH-19001)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NYNRIN (NYN domain and retroviral integrase containing) is a protein-coding gene. GO annotations related to this gene include molecular_function and nucleic acid binding. An important paralog of this gene is ZC3H12C.
Immunogen	amino acids 1476-1604 of Human NYNRIN.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	NYNRIN NYN domain and retroviral integrase containing [Homo sapiens]
Official Symbol	NYNRIN
Synonyms	NYNRIN; NYN domain and retroviral integrase containing; CGIN1; KIAA1305; protein NYNRIN; Cousin of GIN1; protein cousin of GIN1; NYN domain and retroviral integrase catalytic domain-containing protein;
Entrez Gene ID	57523
Protein Refseq	NP_079357.2
UniProt ID	Q9P2P1
Function	molecular_function; nucleic acid binding;