



Anti-HIGD2B (aa 1-47) polyclonal antibody (DPABH-19012)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Multipass membrane protein. Contains 1 HIG1 domain.
Immunogen	antigen sequence: MATLGFVTPE APFESSKPPI FEGLSPTVYS NPEGFKEKFL RKTRENK, corresponding to N terminal amino acids 1-47 of Human HIGD2B
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [HIGD2B HIG1 hypoxia inducible domain family, member 3B \[Homo sapiens \]](#)

Official Symbol	HIGD2B
Synonyms	HIGD2B; HIG1 hypoxia inducible domain family, member 2B; HIGD2BP;
Entrez Gene ID	123346