



Anti-ARMC9 (aa 7-88) polyclonal antibody (DPABH-19034)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ARMC9 (armadillo repeat containing 9) is a protein-coding gene. Diseases associated with ARMC9 include mucopolysaccharidosis iii , and vogt-koyanagi-harada disease.
Immunogen	antigen sequence: HESELLGLVK EYLDFAEFED TLKTFSKECK IKGKPLCKTV GGSFRDSKSL TIQKDLVAAF DNGDQKVFFD LWEEHISSSI RD, corresponding to N terminal amino acids 7-88 of Human ARMC9
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	ARMC9 armadillo repeat containing 10 [Homo sapiens]
Official Symbol	ARMC9
Synonyms	ARMC9; armadillo repeat containing 9; ARM; NS21; KU-MEL-1; lisH domain-containing protein ARMC9; armadillo/beta-catenin-like repeats; melanoma/melanocyte specific protein KU-MEL-1; melanoma/melanocyte-specific tumor antigen KU-MEL-1;
Entrez Gene ID	80210
Protein Refseq	NP_001258395.1
UniProt ID	Q7Z3E5