



# Anti-DLEC1 (aa 186-298) polyclonal antibody (DPABH-19045)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	May act as a tumor suppressor by inhibiting cell proliferation.
Immunogen	antigen sequence, corresponding to internal sequence amino acids 186-298 of Human DLEC1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

Gene Name [DLEC1 deleted in lung and esophageal cancer 2 \[ Homo sapiens \]](#)

<b>Official Symbol</b>	DLEC1
<b>Synonyms</b>	DLEC1; deleted in lung and esophageal cancer 1; F56; DLC1; deleted in lung and esophageal cancer protein 1; DLC-1; deleted in lung and esophageal cancer 1 transcript variant 2;
<b>Entrez Gene ID</b>	<a href="#">9940</a>
<b>Protein Refseq</b>	<a href="#">NP_031361.2</a>
<b>UniProt ID</b>	<a href="#">Q9Y238</a>
<b>Function</b>	molecular_function;