



Anti-BRICD5 (aa 66-169) polyclonal antibody (DPABH-19065)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	BRICD5 (BRICHOS domain containing 5) is a protein-coding gene.
Immunogen	antigen sequence, corresponding to internal sequence amino acids 66-169 of Human C16orf79.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [BRICD5 BRICHOS domain containing 6 \[Homo sapiens \]](#)

Official Symbol	BRICD5
Synonyms	BRICD5; BRICHOS domain containing 5; C16orf79; BRICHOS domain-containing protein 5; BRICHOS domain-containing protein C16orf79;
Entrez Gene ID	283870
Protein Refseq	NP_872369.2
UniProt ID	Q6PL45