



Anti-STARD3NL (aa 1-53) polyclonal antibody (DPABH-19073)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	STARD3NL (STARD3 N-terminal like) is a protein-coding gene. Diseases associated with STARD3NL include parkinsons disease. An important paralog of this gene is STARD3.
Immunogen	antigen sequence: MNHLPEDMEN ALTGSQSSHA SLRNIHSINP TQLMARIESY EGREKKGISD VRR, corresponding to N terminal amino acids 1-53 of Human STARD3NL.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	STARD3NL STARD4 N-terminal like [Homo sapiens]
Official Symbol	STARD3NL
Synonyms	STARD3NL; STARD3 N-terminal like; MENTHO; MLN64 N-terminal domain homolog; STARD3 N-terminal-like protein;
Entrez Gene ID	83930
Protein Refseq	NP_114405.1
UniProt ID	A0A024RA89