



Anti-DEFB129 (aa 53-177) polyclonal antibody (DPABH-19116)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Has antibacterial activity.
Immunogen	antigen, corresponding to amino acids 53-177 of Human DEFB129.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	DEFB129 defensin, beta 130 [Homo sapiens]
-----------	---

Official Symbol	DEFB129
Synonyms	DEFB129; defensin, beta 129; DEFB29; hBD-29; DEFB-29; C20orf87; bA530N10.3; beta-defensin 129; beta defensin 29; beta-defensin 29;
Entrez Gene ID	140881
Protein Refseq	NP_543021.1
UniProt ID	Q9H1M3
Pathway	Beta defensins; Immune System;