



Anti-SNX30 (aa 33-108) polyclonal antibody (DPABH-19147)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in several stages of intracellular trafficking.
Immunogen	antigen sequence: VGGDSTPSPD LLMARSGDK DLILPNGGTP AGTSSPASSS SLLNRLQLDD DIDGETRDLF VIVDDPKKHV CTMETY, corresponding to amino acids 33-108 of Human SNX30.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SNX30 sorting nexin family member 31 [Homo sapiens]
Official Symbol	SNX30
Synonyms	SNX30; sorting nexin family member 30; ATG24A; sorting nexin-30;
Entrez Gene ID	401548
Protein Refseq	NP_001013012.1
UniProt ID	Q5VWJ9
Function	phosphatidylinositol binding; protein binding;