



Anti-C22ORF15 (aa 72-142) polyclonal antibody (DPABH-19179)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C22orf15 (chromosome 22 open reading frame 15) is a protein-coding gene.
Immunogen	antigen sequence: KERAIYVLVR IISKVAQEGE DMASTRYESL LENLDDHYPE LAEELRRLSG LSSVGHNWRK RMGTRRGGRHE QSPTSR, corresponding to amino acids 72-142 of Human C22orf15.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C22ORF15 chromosome 22 open reading frame 16 [Homo sapiens]
Official Symbol	C22ORF15
Synonyms	C22ORF15; chromosome 22 open reading frame 15; N27C7-3; uncharacterized protein C22orf15; protein N27C7-3;
Entrez Gene ID	150248
Protein Refseq	NP_872326.2
UniProt ID	Q8WYQ4