



Rabbit anti-Human C7orf46 polyclonal antibody (DPABH-19221)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FAM221A (family with sequence similarity 221, member A) is a protein-coding gene.
Immunogen	antigen sequence: SGFHSCFTCA CGQPAYAHDT VVETKQERLA QEKVPGQDIP YAAMGGLTGF SSLAEGYMRL DDSGIGVP corresponding to internal sequence amino acids 148-215 of Human C7orf46.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	FAM221A family with sequence similarity 222, member A [Homo sapiens]
Official Symbol	FAM221A
Synonyms	FAM221A; family with sequence similarity 221, member A; C7orf46; protein FAM221A; uncharacterized protein C7orf46;
Entrez Gene ID	340277
Protein Refseq	NP_001120836.1
UniProt ID	A0A024RA57