



Anti-EFCAB10 (aa 67-127) polyclonal antibody (DPABH-19255)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	EFCAB10 (EF-hand calcium binding domain 10) is a protein-coding gene. GO annotations related to this gene include cAMP-dependent protein kinase regulator activity and calcium ion binding.
Immunogen	Recombinant fragment, corresponding to amino acids 67-127 of Human EFCAB10.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	EFCAB10 EF-hand calcium binding domain 11 [Homo sapiens]
Official Symbol	EFCAB10
Synonyms	EFCAB10; EF-hand calcium binding domain 10;
Entrez Gene ID	100130771
Function	calcium ion binding;