



Anti-FAM171A2 (aa 186-281) polyclonal antibody (DPABH-19273)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FAM171A2 (family with sequence similarity 171, member A2) is a protein-coding gene. An important paralog of this gene is FAM171B.
Immunogen	Recombinant fragment corresponding to amino acids 186-281 (RAFPAFLGTE ASSSGNGSWL ELMPLTAVSV HLLTGNNGTEV PLSGPIHLSL PVPSETRALT VGTSIPAWRF DPKSGLWVRN GTGVIRKEGR QLYWTF) of Human FAM171A2 (A8MVW0).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	FAM171A2 family with sequence similarity 171, member A3 [Homo sapiens]
Official Symbol	FAM171A2
Synonyms	FAM171A2; family with sequence similarity 171, member A2; protein FAM171A2;
Entrez Gene ID	284069
Protein Refseq	NP_940877.2
UniProt ID	A8MVW0