



Anti-CCDC171 (aa 542-638) polyclonal antibody (DPABH-19282)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCDC171 (coiled-coil domain containing 171) is a protein-coding gene.
Immunogen	Recombinant fragment corresponding to amino acids 542-638 (KAQAAQSESE LQKLSQAFHK DAEKLTFLH TLYQHLVAGC VLIKQPEGML DKFSWSELCA VLQENVDALI ADLNRANEKI RHLEYICKNK SDTMREL) of Human C9orf93 (Q6TFL3).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CCDC171 coiled-coil domain containing 172 [Homo sapiens]
Official Symbol	CCDC171
Synonyms	CCDC171; coiled-coil domain containing 171; C9orf93; bA536D16.1; bA778P13.1; coiled-coil domain-containing protein 171; myosin tail domain containing protein;
Entrez Gene ID	203238
Protein Refseq	NP_775821.2
UniProt ID	Q6TFL3