



Anti-ARHGAP23 (aa 509-586) polyclonal antibody (DPABH-19303)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.
Immunogen	antigen sequence: MNLGFGDESP EPEASGRGER LGRKVAPLAT TEDSLASIPF IDEPTSPSID LQAKHVPASA VVSSAMNSAP VLGTSPPSS, corresponding to amino acids 509-586 of Human ARHGAP23.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ARHGAP23 Rho GTPase activating protein 24 [Homo sapiens]
Official Symbol	ARHGAP23
Synonyms	ARHGAP23; Rho GTPase activating protein 23; rho GTPase-activating protein 23; rho-type GTPase-activating protein 23;
Entrez Gene ID	57636
Protein Refseq	NP_001186346.1
UniProt ID	Q9P227
Pathway	Rho GTPase cycle; Signaling by Rho GTPases;
Function	GTPase activator activity;