



Anti-PTRHD1 (aa 76-140) polyclonal antibody (DPABH-19313)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PTRHD1 (peptidyl-tRNA hydrolase domain containing 1) is a protein-coding gene. GO annotations related to this gene include aminoacyl-tRNA hydrolase activity and protein tyrosine phosphatase activity.
Immunogen	Recombinant fragment corresponding to amino acids 76-140 (GRMRKVVLEA PDETTLKELA ETLQQKNIDH MLWLEQPENI ATCIALRPYP KEEVGQYLKK FRLFK) of Human PTRHD1 (Q6GMV3).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	PTRHD1 peptidyl-tRNA hydrolase domain containing 2 [Homo sapiens]
Official Symbol	PTRHD1
Synonyms	PTRHD1; peptidyl-tRNA hydrolase domain containing 1; C2orf79; putative peptidyl-tRNA hydrolase PTRHD1; peptidyl-tRNA hydrolase domain-containing protein 1;
Entrez Gene ID	391356
Protein Refseq	NP_001013685.1
UniProt ID	Q6GMV3
Function	aminoacyl-tRNA hydrolase activity; protein tyrosine phosphatase activity;