



Anti-PTCD1 (aa 435-537) polyclonal antibody (DPABH-19336)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PTCD1 (pentatricopeptide repeat domain 1) is a protein-coding gene. GO annotations related to this gene include tRNA binding. An important paralog of this gene is ATP5J2-PTCD1.
Immunogen	antigen, corresponding to amino acids 435-537 of Human PTCD1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [PTCD1 pentatricopeptide repeat domain 2 \[Homo sapiens \]](#)

Official Symbol	PTCD1
Synonyms	PTCD1; pentatricopeptide repeat domain 1; pentatricopeptide repeat-containing protein 1, mitochondrial;
Entrez Gene ID	26024
Protein Refseq	NP_056360.2
UniProt ID	A4D273
Function	poly(A) RNA binding; tRNA binding;