



Anti-C7ORF25 (aa 320-418) polyclonal antibody (DPABH-19347)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C7orf25 (chromosome 7 open reading frame 25) is a protein-coding gene.
Immunogen	antigen sequence: GPGERERATV LIKRINVVPD QPSERALRLV ASSKINSRSL TIFGTGDTLK AITMTANS GF VRAANNQGVK FSVFIHQ PRA LTESKEALAT PLPKDYTTD, corresponding to amino acids 320-418 of Human C7orf25.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C7ORF25 chromosome 7 open reading frame 26 [Homo sapiens]
Official Symbol	C7ORF25
Synonyms	C7ORF25; chromosome 7 open reading frame 25; UPF0415 protein C7orf25;
Entrez Gene ID	79020
Protein Refseq	NP_001093328.1
UniProt ID	Q9BPX7