



Anti-ZC3HAV1L (aa 207-296) polyclonal antibody (DPABH-19379)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ZC3HAV1L (zinc finger CCCH-type, antiviral 1-like) is a protein-coding gene.
Immunogen	antigen sequence: CKRSHQLIHA ASLKLLQDQG LNIPSVVNFQ IISTYKHKML HKMLENTDNS SPSTEHSQGL EKQGVHAAGA AEAGPLASVP AQSAAKKPCPV, corresponding to amino acids 207-296 of Human ZC3HAV1L.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ZC3HAV1L zinc finger CCCH-type, antiviral 2-like [Homo sapiens]
Official Symbol	ZC3HAV1L
Synonyms	ZC3HAV1L; zinc finger CCCH-type, antiviral 1-like; C7orf39; zinc finger CCCH-type antiviral protein 1-like;
Entrez Gene ID	92092
Protein Refseq	NP_542391.2
UniProt ID	Q96H79