



Anti-C7ORF76 (aa 3-63) polyclonal antibody (DPABH-19399)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C7orf76 (chromosome 7 open reading frame 76) is a protein-coding gene.
Immunogen	antigen sequence: CQDSNICAVF AVQGGKVGRK HGIKRGRRPS IRSPAQRARG PWIHESKHPA FAKQQINLEM P, corresponding to N terminal amino acids 3-63 of Human FLJ42280.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C7ORF76 chromosome 7 open reading frame 77 [Homo sapiens]
Official Symbol	C7ORF76
Synonyms	C7ORF76; chromosome 7 open reading frame 76; putative uncharacterized protein C7orf76;
Entrez Gene ID	401388
Protein Refseq	NP_001188379.1
UniProt ID	Q6ZVN7