



Anti-MALSU1 (aa 139-229) polyclonal antibody (DPABH-19438)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	MALSU1 (mitochondrial assembly of ribosomal large subunit 1) is a protein-coding gene. Diseases associated with MALSU1 include pneumonia , and tuberculosis.
Immunogen	antigen sequence: AFYVVKMYKH LKCKRDPHVK IEGKDTDDWL CVDFGSMVIH LMLPETREIY ELEKLWTLRS YDDQLAQIAP ETVPEDFILG IEDDTSSVTP V, corresponding to amino acids 139-229 of Human C7orf30 (Q96EH3).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	MALSU1 mitochondrial assembly of ribosomal large subunit 2 [Homo sapiens]
Official Symbol	MALSU1
Synonyms	MALSU1; mitochondrial assembly of ribosomal large subunit 1; mtRsfA; C7orf30; mitochondrial assembly of ribosomal large subunit protein 1;
Entrez Gene ID	115416
Protein Refseq	NP_612455.1
UniProt ID	Q96EH3