



Anti-C9ORF16 (aa 5-83) polyclonal antibody (DPABH-19476)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	antigen sequence: NGDLGMPVEA GAEGEEDGFG EAEYAAINSM LDQINSCLDH LEEKNDHLHA RLQELLESNR QTRLEFQQQL GEAPSDASP, corresponding to amino acids 5-83 of Human C9orf16.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Concentration	0.27 mg/ml
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol. Note: No stabilizing proteins added.
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C9ORF16 chromosome 9 open reading frame 17 [Homo sapiens]
Official Symbol	C9ORF16
Synonyms	C9ORF16; chromosome 9 open reading frame 16; EST00098; UPF0184 protein C9orf16;
Entrez Gene ID	79095
Protein Refseq	NP_077017.1
UniProt ID	Q9BUW7
Function	molecular_function;