



Anti-C2ORF66 (aa 40-117) polyclonal antibody (DPABH-19509)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C2orf66 (chromosome 2 open reading frame 66) is a protein-coding gene.
Immunogen	antigen sequence, corresponding to amino acids 40-117 (TVRNEDKWKP LNNPRNRDLF FRRLQAYFKG RGLDLGTFPN PFPTNENPRP LSFQSELTAS ASADYEEQKN SFHNYLKG) of Human C2orf66.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	C2ORF66 chromosome 2 open reading frame 67 [Homo sapiens]
Official Symbol	C2ORF66
Synonyms	C2ORF66; chromosome 2 open reading frame 66; UNQ6411; IIDS6411; uncharacterized protein C2orf66;
Entrez Gene ID	401027
Protein Refseq	NP_998773.1
UniProt ID	Q6UXQ4