



# Anti-GPR128 (aa 237-357) polyclonal antibody (DPABH-19598)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                     |                                                                                       |
|---------------------|---------------------------------------------------------------------------------------|
| Antigen Description | Orphan receptor.                                                                      |
| Immunogen           | antigen sequence, corresponding to amino acids 237-357 of Human GPCR GPR128.          |
| Isotype             | IgG                                                                                   |
| Source/Host         | Rabbit                                                                                |
| Species Reactivity  | Human                                                                                 |
| Purification        | Immunogen affinity purified                                                           |
| Conjugate           | Unconjugated                                                                          |
| Applications        | IHC-P                                                                                 |
| Format              | Liquid                                                                                |
| Size                | 100 µl                                                                                |
| Buffer              | pH: 7.20; Constituents: 59% PBS, 40% Glycerol                                         |
| Preservative        | 0.02% Sodium Azide                                                                    |
| Storage             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

## GENE INFORMATION

Gene Name [GPR128 G protein-coupled receptor 129 \[ Homo sapiens \]](#)

|                        |                                                                                  |
|------------------------|----------------------------------------------------------------------------------|
| <b>Official Symbol</b> | GPR128                                                                           |
| <b>Synonyms</b>        | GPR128; G protein-coupled receptor 128; probable G-protein coupled receptor 128; |
| <b>Entrez Gene ID</b>  | <a href="#">84873</a>                                                            |
| <b>Protein Refseq</b>  | <a href="#">NP_116176.2</a>                                                      |
| <b>UniProt ID</b>      | <a href="#">Q6ZMH4</a>                                                           |
| <b>Function</b>        | G-protein coupled receptor activity;                                             |