



Anti-C9ORF172 (aa 307-400) polyclonal antibody (DPABH-19652)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C9orf172 (chromosome 9 open reading frame 172) is a protein-coding gene.
Immunogen	amino acids 307-400 of Human C9orf172.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name	C9ORF172 chromosome 9 open reading frame 173 [Homo sapiens]
-----------	---

Official Symbol	C9ORF172
Synonyms	C9ORF172; chromosome 9 open reading frame 172; uncharacterized protein C9orf172;
Entrez Gene ID	389813
Protein Refseq	NP_001073951.2
UniProt ID	C9J069