



Anti-COX19 (aa 6-80) polyclonal antibody (DPABH-19663)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	COX19 (cytochrome c oxidase assembly homolog 19 (<i>S. cerevisiae</i>)) is a protein-coding gene. Diseases associated with COX19 include menkes disease.
Immunogen	antigen, corresponding to amino acids 6-80 (NFGTKSFQPR PPDKGSFPLD HLGECKSFKE KFMKCLHNNN FENALCRKES KEYLECRMER KLMLQEPLEK LGFGD) of Human COX19.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	COX19 cytochrome c oxidase assembly homolog 20 (S. cerevisiae) [Homo sapiens]
Official Symbol	COX19
Synonyms	COX19; cytochrome c oxidase assembly homolog 19 (S. cerevisiae); cytochrome c oxidase assembly protein COX19; hCOX19; COX19 cytochrome c oxidase assembly homolog;
Entrez Gene ID	90639
Protein Refseq	NP_001026788.1
UniProt ID	Q49B96