



Anti-C7ORF57 (aa 217-292) polyclonal antibody (DPABH-19717)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C7orf57 (chromosome 7 open reading frame 57) is a protein-coding gene. Diseases associated with C7orf57 include lateral sclerosis , and amyotrophic lateral sclerosis.
Immunogen	antigen, corresponding to amino acids 217-292 (ISNGYKDEWL QQQQRADSDK RTPKTSRASV LSQSPRDLEG PQDAARLQDA EASEGPEDTP GPEESVSAST PA) of Human C7orf57.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C7ORF57 chromosome 7 open reading frame 58 [Homo sapiens]
Official Symbol	C7ORF57
Synonyms	C7ORF57; chromosome 7 open reading frame 57; uncharacterized protein C7orf57;
Entrez Gene ID	136288
Protein Refseq	NP_001093629.1
UniProt ID	Q86XB5