



Anti-AAED1 (aa 84-142) polyclonal antibody (DPABH-19728)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	AAED1 (AhpC/TSA antioxidant enzyme domain containing 1) is a protein-coding gene. Diseases associated with AAED1 include gastric cancer.
Immunogen	antigen sequence: SFLQEANVTI IVIGQSSYHH IEPFCKLTGY SHEIYVDPER EIYKRLGMKR GEEIASSGQ, corresponding to amino acids 84-142 of Human C9orf21.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	AAED1 AhpC/TSA antioxidant enzyme domain containing 2 [Homo sapiens]
Official Symbol	AAED1
Synonyms	AAED1; AhpC/TSA antioxidant enzyme domain containing 1; C9orf21; thioredoxin-like protein AAED1; UPF0308 protein C9orf21; ahpC/TSA antioxidant enzyme domain-containing protein 1;
Entrez Gene ID	195827
Protein Refseq	NP_714542.1
UniProt ID	Q7RTV5