



Anti-WDR60 (aa 800-913) polyclonal antibody (DPABH-19739)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RICH1 is a GTPase-activating protein (GAP). GAPs stimulate the intrinsic GTP hydrolysis of small G proteins, such as RHOA (MIM 165390), RAC1 (MIM 602048), and CDC42 (MIM 116952).[supplied by OMIM, Apr 2004]
Immunogen	antigen sequence, corresponding to amino acids 800-913 of Human WDR60.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	WDR60 WD repeat domain 61 [Homo sapiens]
Official Symbol	WDR60
Synonyms	WDR60; WD repeat domain 60; SRPS6; SRTD8; FAP163; WD repeat-containing protein 60;
Entrez Gene ID	55112
Protein Refseq	NP_060521.4
UniProt ID	Q8WVS4