



# Anti-FAM71F1 (aa 200-268) polyclonal antibody (DPABH-19824)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	FAM71F1 (family with sequence similarity 71, member F1) is a protein-coding gene. Diseases associated with FAM71F1 include obesity. An important paralog of this gene is FAM71E1.
<b>Immunogen</b>	antigen sequence: GLSITTKDPR ILVTHCLVPK NCSSPSGDSK LVQKKLQASQ PSESLIQLMT KGESEALSQI FADLHQQNQ, corresponding to amino acids 200-268 of Human FAM71F1.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P, WB
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">FAM71F1 family with sequence similarity 71, member F2 [ Homo sapiens ]</a>
<b>Official Symbol</b>	FAM71F1
<b>Synonyms</b>	FAM71F1; family with sequence similarity 71, member F1; FAM137A; NYD-SP18; protein FAM71F1; testes development-related NYD-SP18; testis development protein NYD-SP18; family with sequence similarity 137, member A;
<b>Entrez Gene ID</b>	<a href="#">84691</a>
<b>Protein Refseq</b>	<a href="#">NP_001269717.1</a>
<b>UniProt ID</b>	<a href="#">Q96KD3</a>