



Rabbit anti-Human C9orf116 polyclonal antibody (DPABH-19854)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C9orf116 (chromosome 9 open reading frame 116) is a protein-coding gene.
Immunogen	antigen sequence: PTVHEMPKVF YPNSNKFSQQ LAAGGMFRNN TLNVYLEKSI VTGPDNCITS CDRLNFHPSY NINRPSICD, corresponding to C terminal amino acids 68-136 of Human C9orf116.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C9ORF116 chromosome 9 open reading frame 117 [Homo sapiens]
Official Symbol	C9ORF116
Synonyms	C9ORF116; chromosome 9 open reading frame 116; PIERCE1; RbEST47; UPF0691 protein C9orf116; pierce 1; p53-induced expression in RB-null cells protein 1;
Entrez Gene ID	138162
Protein Refseq	NP_001041730.1
UniProt ID	Q5BN46