



Anti-SVEP1 (aa 736-827) polyclonal antibody (DPABH-19887)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role in the cell attachment process.
Immunogen	antigen sequence: GDFICTPDNT GVNCTLTCLG GYDFTEGSTD KYCAYEDGV WKPTYTTEWP DCAKKRFANH GFKSFEMFYK AARCDDTDLM KKFSEAFETT LG, corresponding to internal sequence amino acids 736-827 of Human SVEP1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SVEP1 sushi, von Willebrand factor type A, EGF and pentraxin domain containing 2 [Homo sapiens]
Official Symbol	SVEP1
Synonyms	SVEP1; sushi, von Willebrand factor type A, EGF and pentraxin domain containing 1; CCP22; SELOB; SEL-OB; C9orf13; POLYDOM; sushi, von Willebrand factor type A, EGF and pentraxin domain-containing protein 1; selectin-like protein; CCP module-containing protein 22; selectin-like osteoblast-derived protein; serologically defined breast cancer antigen NY-BR-38;
Entrez Gene ID	79987
Protein Refseq	NP_699197.3
UniProt ID	B3KQM1
Function	calcium ion binding; chromatin binding;