



# Anti-ANKFN1 (aa 336-427) polyclonal antibody (DPABH-19948)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                                                                                  |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antigen Description</b> | ANKFN1 (ankyrin-repeat and fibronectin type III domain containing 1) is a protein-coding gene. Diseases associated with ANKFN1 include cannabis dependence.                                      |
| <b>Immunogen</b>           | antigen sequence: GLTMGQQYFV QVSAYNMKGW GPAQTTTPAC ASPSNWKDYD DREPRHKGQS EVLEGLLQQV RALHQHYSCR ESTKLQTTGR KQSVSRSLKH LF, corresponding to internal sequence amino acids 336-427 of Human ANKFN1. |
| <b>Isotype</b>             | IgG                                                                                                                                                                                              |
| <b>Source/Host</b>         | Rabbit                                                                                                                                                                                           |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                                            |
| <b>Purification</b>        | Immunogen affinity purified                                                                                                                                                                      |
| <b>Conjugate</b>           | Unconjugated                                                                                                                                                                                     |
| <b>Applications</b>        | IHC-P, ICC/IF                                                                                                                                                                                    |
| <b>Format</b>              | Liquid                                                                                                                                                                                           |
| <b>Size</b>                | 100 µl                                                                                                                                                                                           |
| <b>Buffer</b>              | pH: 7.20; Constituents: 59% PBS, 40% Glycerol                                                                                                                                                    |
| <b>Preservative</b>        | 0.02% Sodium Azide                                                                                                                                                                               |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.                                                                                                            |

## GENE INFORMATION

|                        |                                                                                                                                           |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>       | <a href="#">ANKFN1 ankyrin-repeat and fibronectin type III domain containing 2 [ Homo sapiens ]</a>                                       |
| <b>Official Symbol</b> | ANKFN1                                                                                                                                    |
| <b>Synonyms</b>        | ANKFN1; ankyrin-repeat and fibronectin type III domain containing 1; ankyrin repeat and fibronectin type-III domain-containing protein 1; |
| <b>Entrez Gene ID</b>  | <a href="#">162282</a>                                                                                                                    |
| <b>Protein Refseq</b>  | <a href="#">NP_694960.2</a>                                                                                                               |
| <b>UniProt ID</b>      | <a href="#">Q8N957</a>                                                                                                                    |