



Rabbit anti-Human FAM205A polyclonal antibody (DPABH-19981)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This locus represents a transcribed pseudogene of a nearby locus on chromosome 7, which encodes a putative methyltransferase. There is also a third closely related pseudogene locus in this region. There is extensive alternative splicing at this locus.
Immunogen	antigen sequence, corresponding to internal sequence amino acids 922-1031 of Human FAM205A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	FAM205A family with sequence similarity 206, member A [Homo sapiens]
Official Symbol	FAM205A
Synonyms	FAM205A; family with sequence similarity 205, member A; C9orf144B; protein FAM205A; transmembrane protein C9orf144B;
Entrez Gene ID	259308
Protein Refseq	NP_001135389.1
UniProt ID	B4DY95