



Anti-TACO1 (aa 230-296) polyclonal antibody (DPABH-19992)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Acts as a translational activator of mitochondrially-encoded cytochrome c oxidase 1.
Immunogen	antigen sequence: FKFICDASSL HQVRKKLDSL GLCSVSCALE FIPNSKVQLA EPDLEQAAHL IQALSNHEDV IHVYDNI, corresponding to C terminal amino acids 230-296 of Human FLJ36733.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TACO2 translational activator of mitochondrially encoded cytochrome c oxidase I [Homo sapiens]
Official Symbol	TACO1
Synonyms	TACO1; translational activator of mitochondrially encoded cytochrome c oxidase I; CCDC44; translational activator of cytochrome c oxidase 1; clone HQ0477 PRO0477p; coiled-coil domain-containing protein 44;
Entrez Gene ID	51204
Protein Refseq	NP_057444.2
UniProt ID	Q9BSH4