



Anti-INTS1 (aa 289-379) polyclonal antibody (DPABH-20010)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Component of the Integrator complex, a complex involved in the small nuclear RNAs (snRNA) U1 and U2 transcription and in their 3-box-dependent processing. The Integrator complex is associated with the C-terminal domain (CTD) of RNA polymerase II largest subunit (POLR2A) and is recruited to the U1 and U2 snRNAs genes.
Immunogen	antigen sequence: LTEEEDSQTE LLIAEEKLSP EQEQLMPRY EELAESVEEY VLDMLRDQLN RRQPIDNVSRL LLRLLTSTC GYKEVRLAV QKLEMWLQNP K, corresponding to amino acids 289-379 of Human INTS1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	INTS1 integrator complex subunit 2 [Homo sapiens]
Official Symbol	INTS1
Synonyms	INTS1; integrator complex subunit 1; INT1; NET28;
Entrez Gene ID	26173
Protein Refseq	NP_001073922.2
UniProt ID	Q8N201