



Anti-UBAP2 (aa 429-533) polyclonal antibody (DPABH-20011)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a phosphoprotein that interacts with several members of the E2F family of proteins. The protein localizes to the nucleus, and is present throughout the cell cycle except during mitosis. It functions to modulate E2F-regulated transcription and stimulate proliferation.
Immunogen	antigen, corresponding to amino acids 429-533 of Human UBAP2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	UBAP2 ubiquitin associated protein 3 [Homo sapiens]
Official Symbol	UBAP2
Synonyms	UBAP2; ubiquitin associated protein 2; UBAP-2; ubiquitin-associated protein 2; AD-012 protein;
Entrez Gene ID	55833
Protein Refseq	NP_001269458.1
UniProt ID	B4EQR5
Function	poly(A) RNA binding;