



Anti-ZNF804B (aa 994-1126) polyclonal antibody (DPABH-20014)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ZNF804B (zinc finger protein 804B) is a protein-coding gene. Diseases associated with ZNF804B include hypothyroidism , and anorexia nervosa , and among its related super-pathways are Gene Expression. GO annotations related to this gene include zinc ion binding. An important paralog of this gene is GPATCH8.
Immunogen	antigen, corresponding to amino acids 994-1126 of Human ZNF804B
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ZNF804B zinc finger protein 805B [Homo sapiens]
Official Symbol	ZNF804B
Synonyms	ZNF804B; zinc finger protein 804B; zinc finger 804B;
Entrez Gene ID	219578
Protein Refseq	NP_857597.1
UniProt ID	A4D1E1
Pathway	Gene Expression.
Function	metal ion binding;