



Anti-LRRC46 (aa 231-315) polyclonal antibody (DPABH-20017)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	LRRC46 (leucine rich repeat containing 46) is a protein-coding gene. An important paralog of this gene is LRRC6.
Immunogen	antigen sequence: PTLTDLPLLP GVPMAGDSSP SATPAQGEET VPEAVSSPQA SSPTKKPCSL IPRGHQSSFW GRKGARAATA PKASVAEAPS TTKTT, corresponding to amino acids 231-315 of Human LRRC46
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	LRRC46 leucine rich repeat containing 47 [Homo sapiens]
Official Symbol	LRRC46
Synonyms	LRRC46; leucine rich repeat containing 46; leucine-rich repeat-containing protein 46;
Entrez Gene ID	90506
Protein Refseq	NP_219481.1
UniProt ID	Q96FV0