



Anti-CYSRT1 (aa 77-144) polyclonal antibody (DPABH-20020)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein Epitope Signature Tag: PIQNQQAWQQ PGNPYSSSQ R QAGLTYAGPP PVGRGDDIAH HCCCCPCCHC CHCPPFCRCH SCCCCVIS, corresponding to amino acids 77 - 144 of Human C9orf169.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [CYSRT1 cysteine-rich tail protein 2 \[Homo sapiens \]](#)

Official Symbol	CYSRT1
Synonyms	CYSRT1; cysteine-rich tail protein 1; C9orf169; UPF0574 protein C9orf169;
Entrez Gene ID	375791
Protein Refseq	NP_945352.3
UniProt ID	A8MQ03