



Anti-QRIC2 (aa 1450-1554) polyclonal antibody (DPABH-20022)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	QRICH2 (glutamine rich 2) is a protein-coding gene.
Immunogen	antigen sequence, corresponding to amino acids 1450-1554 of Human QRICH2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [QRICH2 glutamine rich 3 \[Homo sapiens \]](#)

Official Symbol	QRICH2
Synonyms	QRICH2; glutamine rich 2; glutamine-rich protein 2;
Entrez Gene ID	84074
Protein Refseq	NP_115510.1
UniProt ID	Q9H0J4