



Anti-CDR2L (aa 395-464) polyclonal antibody (DPABH-20025)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CDR2L (cerebellar degeneration-related protein 2-like) is a protein-coding gene. Diseases associated with CDR2L include cerebellar degeneration. An important paralog of this gene is CDR2.
Immunogen	antigen sequence corresponding to amino acids 395-464 (DSSWRDLRGG EEGQGEVKAG EKSLSQHVEA VDKRLEQSQP EYKALFKEIF SRIQKTKADI NATKVKTHSS) of Human CDRL2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CDR2L cerebellar degeneration-related protein 3-like [Homo sapiens]
Official Symbol	CDR2L
Synonyms	CDR2L; cerebellar degeneration-related protein 2-like; HUMPPA; paraneoplastic antigen; paraneoplastic 62 kDa antigen;
Entrez Gene ID	30850
Protein Refseq	NP_055418.2
UniProt ID	A0A024R8R6