



Anti-SPATA20 (aa 52-139) polyclonal antibody (DPABH-20026)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role in fertility regulation.
Immunogen	antigen sequence: SHPSSTPQRV PNRLIHEKSP YLLQHAYNPV DWYPWGQEAF DKARKENKPI FLSVGYSTCH WCHMMEEESF QNEEIGRLLS EDFVSVKV, corresponding to amino acids 52-139 of Human SPATA20.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SPATA20 spermatogenesis associated 21 [Homo sapiens]
Official Symbol	SPATA20
Synonyms	SPATA20; spermatogenesis associated 20; SSP411; Tisp78; spermatogenesis-associated protein 20; sperm protein SSP411; sperm-specific protein 411; transcript increased in spermiogenesis 78;
Entrez Gene ID	64847
Protein Refseq	NP_001245301.1
UniProt ID	Q8TB22
Function	catalytic activity;