



Anti-NKTR (aa 1061-1168) polyclonal antibody (DPABH-20027)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Component of a putative tumor-recognition complex. Involved in the function of NK cells.
Immunogen	antigen sequence corresponding to amino acids 1061-1168 of Human NK-TR protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [NKTR natural killer cell triggering receptor \[Homo sapiens \]](#)

Official Symbol	NKTR
Synonyms	NKTR; natural killer cell triggering receptor; p104; NK-tumor recognition protein; NK-TR protein; natural killer triggering receptor; natural killer-tumor recognition sequence; natural-killer cells cyclophilin-related protein;
Entrez Gene ID	4820
Protein Refseq	NP_005376.2
UniProt ID	A0A024R2N6
Function	cyclosporin A binding; peptidyl-prolyl cis-trans isomerase activity;