



Anti-KIF27 (aa 553-668) polyclonal antibody (DPABH-20028)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Plays an essential role in motile ciliogenesis.
Immunogen	antigen sequence, corresponding to amino acids 553-668 of Human KIF27.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [KIF27 kinesin family member 28 \[Homo sapiens \]](#)

Official Symbol	KIF27
Synonyms	KIF27; kinesin family member 27; kinesin-like protein KIF27;
Entrez Gene ID	55582
Protein Refseq	NP_001258856.1
UniProt ID	Q86VH2
Function	ATP binding; microtubule binding; microtubule motor activity;