



Anti-NYAP1 (aa 549-648) polyclonal antibody (DPABH-20030)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NYAP1 (neuronal tyrosine-phosphorylated phosphoinositide-3-kinase adaptor 1) is a protein-coding gene. Diseases associated with NYAP1 include neuronitis. An important paralog of this gene is NYAP2.
Immunogen	antigen sequence, corresponding to amino acids 549-648 (WTYPATAAGL KRPPAYESLK AGGVNLKGC G VGAPSPMVKI QLQEQGTDGG AFASISCAHV IASAGTPEEE EEEVGAATFG AGWALQRKVL YGGRKAKELD) of Human C7orf51.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	NYAP1 neuronal tyrosine-phosphorylated phosphoinositide-3-kinase adaptor 2 [Homo sapiens]
Official Symbol	NYAP1
Synonyms	NYAP1; neuronal tyrosine-phosphorylated phosphoinositide-3-kinase adaptor 1; C7orf51; neuronal tyrosine-phosphorylated phosphoinositide-3-kinase adapter 1;
Entrez Gene ID	222950
Protein Refseq	NP_775835.2
UniProt ID	Q6ZVC0