



Anti-GARNL3 (aa 503-598) polyclonal antibody (DPABH-20031)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	GARNL3 (GTPase activating Rap/RanGAP domain-like 3) is a protein-coding gene. Diseases associated with GARNL3 include bipolar disorder. GO annotations related to this gene include small GTPase regulator activity and GTPase activator activity. An important paralog of this gene is SIPA1.
Immunogen	antigen sequence, corresponding to amino acids 503-598 (LVSTDAGVLL VDDDLPSVPV FDRTPVKQM HVLETLDLLV LRADKGGK DAR LFVFRSLAQ KGLEGKQAGK SRSDCRENKL ETKGCHLYA INTTHS) of Human GARNL3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	GARNL3 GTPase activating Rap/RanGAP domain-like 4 [Homo sapiens]
Official Symbol	GARNL3
Synonyms	GARNL3; GTPase activating Rap/RanGAP domain-like 3; bA356B19.1; GTPase-activating Rap/Ran-GAP domain-like protein 3; GTPase activating RANGAP domain-like 3;
Entrez Gene ID	84253
Protein Refseq	NP_001273708.1
UniProt ID	Q5VWW2
Function	GTPase activator activity; small GTPase regulator activity;