



Anti-ANKFN1 (aa 659-748) polyclonal antibody (DPABH-20034)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ANKFN1 (ankyrin-repeat and fibronectin type III domain containing 1) is a protein-coding gene. Diseases associated with ANKFN1 include cannabis dependence.
Immunogen	antigen sequence, corresponding to amino acids 659-748 (LSGSESMESV DHTSDCPMQL FFYELQMAVK ALLQQINIPL HQARNFRLYT QEVLEMGHNV SFLLLLPASD DVCTAPGQNN PYTPHSGFLN) of Human ANKFN1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ANKFN1 ankyrin-repeat and fibronectin type III domain containing 2 [Homo sapiens]
Official Symbol	ANKFN1
Synonyms	ANKFN1; ankyrin-repeat and fibronectin type III domain containing 1; ankyrin repeat and fibronectin type-III domain-containing protein 1;
Entrez Gene ID	162282
Protein Refseq	NP_694960.2
UniProt ID	Q8N957