



# Anti-FAM187A (aa 19-112) polyclonal antibody (DPABH-20036)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	FAM187A (family with sequence similarity 187, member A) is a protein-coding gene. An important paralog of this gene is FAM187B.
<b>Immunogen</b>	antigen sequence, corresponding to amino acids 19-112 (EIVEKENIFQ RTPCPAFLMF ENAAYLADMS FELPCHCKPE EVPAVVWFYQ KHLGSSHTKV LTDFDGRVLT EAAQVRVGSD MLTRFSIRMF SLLV) of Human FAM187A.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ICC/IF, IHC-P
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">FAM187A family with sequence similarity 188, member A [ Homo sapiens ]</a>
<b>Official Symbol</b>	FAM187A
<b>Synonyms</b>	FAM187A; family with sequence similarity 187, member A; RGD1307851; Ig-like V-type domain-containing protein FAM187A; ig-like V-type domain-containing protein ENSP00000329499 homolog;
<b>Entrez Gene ID</b>	<a href="#">287741</a>
<b>Protein Refseq</b>	<a href="#">NP_001013883.1</a>
<b>UniProt ID</b>	<a href="#">Q6AXV7</a>