



Anti-KRBA2 (aa 320-403) polyclonal antibody (DPABH-20041)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KRBA2 (KRAB-A domain containing 2) is a protein-coding gene. GO annotations related to this gene include nucleic acid binding. An important paralog of this gene is ZBED5.
Immunogen	Recombinant Protein Epitope Signature Tag: SVLDSDSGVE FTNQVVHELN ELWPDLKIVS GKYHPGQSQG SLEGASRDVK NMISTWMQSN HSCHWAKGLR FMQMVRNQAF DVSL, corresponding to amino acids 320 - 403 of Human KRBA2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	KRBA2 KRAB-A domain containing 3 [Homo sapiens]
Official Symbol	KRBA2
Synonyms	KRBA2; KRAB-A domain containing 2; KRAB-A domain-containing protein 2; KRAB A domain containing 2;
Entrez Gene ID	124751
Protein Refseq	NP_998762.1
UniProt ID	Q6ZNG9
Function	nucleic acid binding;